

WPC 80 Bioactive Protein Comparison

<u>LCMS Results</u>	<u>UV pasteurized</u>	<u>A</u>	<u>B</u>
Immunoglobulin	204,450,000	93,655,000	55,602,000
Lactoferrin	108,350,000	14,466,000	6,789,400
Lactoperoxidase	24,395,000	16,863,000	2,483,600
Polymeric immunoglobulin receptor	28,575,000	3,135,400	978,610

Multiple of Others

Immunoglobulin	2.2	3.7
Lactoferrin	7.5	16.0
Lactoperoxidase	1.4	9.8
Polymeric immunoglobulin receptor	9.1	29.2

